




The
University
Of
Sheffield.

Ultraviolet / Visible Spectroscopy Service Sample Submission Form

Department of Chemistry
Dainton Building
Sheffield, UK, S7 3HF



User Information

Company Name:

Address:

.....

Contact:

Email:

Office Phone:

Purchase order No:

Date:

Properties & Hazards

Storage:

Bench Fridge Freezer

Other stability considerations if any:

.....

Hazards:

COSHH No

Toxic Harmful Corrosive

Irritant Explosive Flammable

Other

.....

PLEASE INCLUDE MSDS IF AVAILABLE

Sample Information

- Sample Reference No(s):
- Approx. / exact* sample weight (* delete as appropriate):
.....
- Scan range (Note: 190nm – 1100nm is the full range but can be limited depending on solvent):
.....
- Solvent (Note: for the best spectra in the lower wavelength range choose a solvent with a low UV cut-off value such as water, acetonitrile, methanol or hexane):
- Should the exact sample concentration be reported? Yes No
- Samples will be diluted until all peaks are on scale unless otherwise instructed here:
.....

Structure(s)

Service use only :

Instrument

Filename

Total time