




The  
University  
Of  
Sheffield.

## Ultra-Violet / Visible Spectroscopy Service Sample Submission Form



### User Information

Name: .....

Department: .....

Group: .....

Lab: .....

Email: .....

Office Phone: .....

Charge Code: .....

Date: .....

### Properties & Hazards

#### Storage:

Bench  Fridge  Freezer

#### Other stability considerations if any:

.....

#### Hazards:

COSHH No .....

Toxic  Harmful  Corrosive

Irritant  Explosive  Flammable

Other

.....

### Sample Information

- Sample Reference No(s): .....
- Approx. / exact\* sample weight (\* delete as appropriate):  
.....
- Scan range (Note:190nm – 1100nm is the full range but can be limited depending on solvent):  
.....
- Solvent (Note: for the best spectra in the lower wavelength range choose a solvent with a low UV cut-off value such as water, acetonitrile, methanol or hexane):.....
- Should the exact sample concentration be reported?  Yes  No
- Samples will be diluted until all peaks are on scale unless otherwise instructed here:  
.....

#### Structure(s)

- Return un-used sample ?  Yes  No

Service use only :  
Instrument  
Filename

Total time